

## Station A: MagMAX Viral/Pathogen Nucleic Acid Isolation (v2)


### Code parameters:

- Change the sample number on line 14 (default is 8, max is 96)
- Change the sample volume (µL) on line 15 (default is 400µL)
- Tip rack tracking can be changed to from False to True on line 16 (default is False)

### Pipettes:

- P1000 single-channel on the right mount
- P20 multi-channel on the left mount

### Deck Layout:



**Labware and module requirements:**

- $\leq 96 \times 1.5$  - 2mL tubes [**input - samples**]
- 4 x 24-Position Tube Rack [**holds 1.5/2mL sample tubes**]
- 3 x 1000 $\mu$ L Filter Tip Racks
- 1 x 6-Position Tube Rack [**holds 50mL tube with Binding Bead Mix**]
- 2 x 50mL Falcon Tube [**1 tube per 48 samples**]
- 1 x 2mL Deep Well Plate [**output - destination plate**]
- 1 x 20 $\mu$ L Filter Tip Racks
- 1 x 96-Well Aluminum Block [**holds 1 PCR strip tube in column 1**]
- 1 x 200 $\mu$ L PCR strip tube [**holds Proteinase K**]

**Volume requirements:**

**Note** ~ the below volumes account for a dead volume; the dead volume can be adjusted depending on the calibration of the pipette to the labware, but we've found it's best to have a dead volume of at least 10%


Reagent	Volume per sample	Volume for 8 samples	Volume for 48 samples	Volume for 96 samples
<b>Binding Bead Mix</b>	550 $\mu$ L	4.84mL	29mL	58.1mL
<b>Proteinase K</b>	10 $\mu$ L	90 $\mu$ L	530 $\mu$ L	1060 $\mu$ L

Before you begin:


1. Load the samples onto the tube racks, starting with slot 2, then slot 3, slot 5, and slot 6
2. Create the **Binding Bead Mix for slot 7**

Reagent	Volume per sample	Volume for 8 samples	Volume for 48 samples	Volume for 96 samples
<b>Binding Solution</b>	530 $\mu$ L	4.66mL	28mL	56mL
<b>Total Nucleic Acid Magnetic Beads</b>	20 $\mu$ L	176 $\mu$ L	1050 $\mu$ L	2100 $\mu$ L

3. Add the binding bead mix in the 50mL tube(s) to the rack in **Slot 7**



4. Divide the total volume of the Proteinase K mix by 8 and add it to a 200µL PCR strip tube to be placed in **Slot 1** on top of a 96-well Aluminum Block in column 1



The final destination will be the deep well plate in slot 4. Once the run is finished, move the deep well plate to Station B to complete the remainder of the extraction protocol.